

STUDY OF MODEL LIPID BIOMEMBRANES WITH HIGH-RESOLUTION ADIABATIC SCANNING CALORIMETRY

INTRODUCTION

Scanning calorimetric techniques allow to obtain continuously the evolution of the heat capacity at constant pressure $C_p(T)$ in terms of the power P and the rate dT/dt .

$$C_p = \frac{dQ}{dT} = \frac{dQ/dt}{dT/dt} = \frac{P}{\dot{T}}$$

In an Adiabatic Scanning Calorimeter (ASC) a constant power P is supplied to the sample and the resulting change in temperature $T(t)$ is measured as a function of time from which the rate dT/dt can be calculated. Combining the rate with the constant power results in $C_p(T)$.

$$H(T) = P[t(T) - t(T_0)]$$


Moreover, the enthalpy $H(T)$ is easily obtained from the product of the power P and the time laps between the start of the run at t_0 and the time at which $T(t)$ was reached.

SAMPLES USED

Lipids are the main constituents of cell membranes. Many of these lipids exhibit phase transitions near typical body temperatures, making these transitions biologically relevant.

Here, dilute solutions of a model lipid, 1,2-dipalmitoyl-sn-glycerol-3-phosphocholine (DPPC), were measured at increasing dilution level. Four concentrations were prepared by diluting the same initial stock solution. From the respective solutions, an amount was placed in a commercial high-pressure DSC cell. A Peltier-based ASC was used, and an average heating rate of 1 K/min was achieved. The following samples were studied:

Concentration DPPC	Amount of solution in the cell	Effective amount of lipid in the cell
1%	63.2 mg	632 μg
0.25%	63.7 mg	159 μg
0.08%	71.6 mg	57 μg
0.03%	65.4 mg	20 μg





(a) Concentration dependence of C_p profile. (b) Enthalpy of the ripple to liquid crystalline phase transition in sample with 1% DPPC concentration.

RESULTS

Two phase transitions were detected for all concentrations, although the smaller transition becomes very faint at the lowest concentrations. The phase transitions separate, with increasing temperature, the gel phase L_{β}' , the ripple phase P_{β}' and the liquid crystalline phase L_{α} [1].


The results are better than those [1] of a Privalov-type DSC calorimeter [2], the standard calorimeter type used in the study of solutions of biological systems. In this case, a sample similar to our second concentration (with 159 μg DPPC) was measured, but with a total sample of 500 mg and effective mass of DPPC of 395 μg . ASC shows the transition to be sharper and more symmetric than the DSC.



C_p profile. (a) This work. (b) Ref. 1.

REFERENCES

- [1] Pfeiffer, H., Klose, G., Heremans, K., & Glorieux, C. (2006), Thermotropic phase behaviour of the pseudobinary mixtures of DPPC/ $C_{12}E_5$ and DMPC/ $C_{12}E_5$ determined by differential scanning calorimetry and ultrasonic velocimetry. *Chemistry and Physics of Lipids*, 139(1), 54–67.
- [2] Privalov, P. L. (1980). Scanning microcalorimeters for studying macromolecules. *Pure and Applied Chemistry*, 52(2), 479–497.



Detailed C_p profile. of the ripple to liquid crystalline phase transition. (a) This work. (b) Ref. 1.